

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant

O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Disinfectant | | | | | | |
|--------------------|---|---------------------|---------------------|----------------------|----------------------|----------------|---------------------|------------------------|
| | | Cidex | Cidex OPA | Gigasept AF | Gigasept FF | Virkon | Cidex Plus | Wavicide-01 |
| Manufacturer | | Johnson and Johnson | Johnson and Johnson | Schülke&Mayr UK Ltd. | Schülke&Mayr UK Ltd. | Dupont (Antec) | Johnson and Johnson | Medical Chemical Corp. |
| Level | | HLD | HLD | - | HLD | HLD | HLD | HLD |
| Active Ingredients | | Glutaraldehyde | Orthophthalaldehyde | Succinc dialdehyde | Succinc dialdehyde | Potassium | Glutaraldehyde | Glutaraldehyde |
| FDA Approval | | Y | Y | Y | Y | Y | Y | Y |
| Disinfectant Type | | Liquid | Liquid | Liquid | Liquid | Powder | Liquid | Liquid |
| Convex | C1-6/C1-6T/C1-6i/SC1-6 /SC1-4HS/SC1-4H/ SC1-6H/ SC2-9H/SC1-7H | Δ | O | Δ | X | O | O | O |
| | C1-6CT | Δ | O | Δ | X | Δ | Δ | O |
| Micro-Convex | C5-8 | O | O | O | O | Δ | O | O |
| | C5-8N/C5-8NT/SC2-11H | Δ | O | Δ | X | O | O | O |
| Linear | L3-8/L3-8i/L3-12/L3-12T /L3-12i/L3-12H/ L3-12H ^{WD} /L3-2X/L8-17 /L8-17i/L8-17X/L3-15H | Δ | O | Δ | X | O | O | O |
| | L10-25H/SL3-19H /SL3-19X | Δ | O | Δ | O | O | O | Δ |
| Phased | MP1-5X/ SP1-5/SP1-5i/ SP1-5X/P1-5CT/SP1-5T/ SP3-8/SP3-8T | Δ | O | Δ | X | O | O | O |
| Volume | VC1-6/VC1-6T/SVC1-6/ SVC1-6H/VE3-10H/ SVC1-8H | Δ | O | Δ | O | O | O | Δ |
| Endo | EV3-10/EV3-10T/ EV3-10H/EV3-10X | Δ | O | Δ | X | O | O | O |
| | EC3-10/EC3-10T/ EC3-10H/EC3-10X/ EC(V)2-11H | Δ | O | Δ | X | O | O | O |
| | EN3-10 | O | O | Δ | Δ | Δ | O | O |
| | E3-10 | Δ | O | Δ | Δ | O | O | Δ |
| | E3-10H | Δ | O | Δ | Δ | O | O | Δ |
| IO | IO3-12/IO8-17/IO8-17T | Δ | O | Δ | X | O | O | O |
| CW | CW2.0/CW5.0 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| TEE | TEE3-7 | Δ | O | Δ | Δ | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | Δ | O | Δ | X | O | O | O |
| Others | US-Guider | Δ | O | Δ | X | X | O | O |
| | Strain relief&Cable | Δ | O | Δ | X | O | O | X |

※ CAUTION ※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

***NOTE**

- Antigermix, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant

O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Disinfectant | | | | | | |
|--------------------|---|------------------|---------------------|------------------------|----------------|------------|----------------|----------------|
| | | AIDAL PLUS | Ster-Bac Blu | Hibitane | Metricide 14 | Sporicidin | Omnicide -FG2 | Nuclear |
| Manufacturer | | Whiteley Medical | Ecolab | COOPERS | Metrex | CONTEC | Interhatch | Ecolab |
| Level | | HLD | - | - | HLD | HLD | - | - |
| Active Ingredients | | Glutaraldehyde | Quaternary Ammonium | Clorhexidine gluconate | Glularaldehyde | Phenol | Glutaraldehyde | Glutaraldehyde |
| FDA Approval | | Y | N | N | Y | Y | Y | N |
| Disinfectant Type | | Liquid | Liquid | Liquid | Liquid | Liquid | Liquid | Liquid |
| Convex | C1-6/C1-6T/C1-6i/SC1-6 /SC1-4HS/SC1-4H/ SC1-6H/ SC2-9H/SC1-7H | O | Δ | Δ | Δ | X | Δ | Δ |
| | C1-6CT | O | Δ | Δ | Δ | O | Δ | Δ |
| Micro-Convex | C5-8 | Δ | O | O | O | Δ | Δ | Δ |
| | C5-8N/C5-8NT/SC2-11H | O | Δ | Δ | Δ | X | Δ | Δ |
| Linear | L3-8/L3-8i/L3-12/L3-12T /L3-12i/L3-12H/ L3-12H ^{WD} /L3-2X/L8-17 /L8-17i/L8-17X/L3-15H | O | Δ | Δ | Δ | X | Δ | Δ |
| | L10-25H/SL3-19H /SL3-19X | O | Δ | Δ | Δ | Δ | Δ | Δ |
| Phased | MP1-5X/ SP1-5/SP1-5i/ SP1-5X/P1-5CT/SP1-5T/ SP3-8/SP3-8T | O | Δ | Δ | Δ | X | Δ | Δ |
| Volume | VC1-6/VC1-6T/SVC1-6/ SVC1-6H/VE3-10H/ SVC1-8H | O | Δ | Δ | Δ | Δ | Δ | Δ |
| Endo | EV3-10/EV3-10T/ EV3-10H/EV3-10X | O | Δ | Δ | Δ | X | Δ | Δ |
| | EC3-10/EC3-10T/ EC3-10H/EC3-10X/ EC(V)2-11H | O | Δ | Δ | Δ | X | Δ | Δ |
| | EN3-10 | Δ | Δ | Δ | Δ | Δ | O | O |
| | E3-10 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| | E3-10H | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| IO | IO3-12/IO8-17/IO8-17T | O | Δ | Δ | Δ | X | Δ | Δ |
| CW | CW2.0/CW5.0 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| TEE | TEE3-7 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | O | Δ | Δ | Δ | X | Δ | Δ |
| Others | US-Guider | O | Δ | Δ | Δ | O | Δ | Δ |
| | Strain relief&Cable | X | Δ | Δ | Δ | X | Δ | Δ |

※ CAUTION ※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

***NOTE**

- Antigermix, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant

O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Disinfectant | | | | | | |
|--------------------|---|---------------------|-----------|------------------|---------------------|-------------------------|---------------------|------------------------|
| | | Sani-Cloth HB | PeraSafe | Super Sani-Cloth | Cleanisept wipes | CLEANISEPT® WIPES forte | Protex | Tristel Duo |
| Manufacturer | | NICE PAK | VWR | NICE PAK | Dr. Schumacher GmbH | Dr. Schumacher GmbH | PARKER | Tristel Solutions Ltd. |
| Level | | LLD | HLD | - | LLD | LLD | LLD | HLD |
| Active Ingredients | | Quaternary Ammonium | Peracetic | Alcohol | Quaternary Ammonium | Quaternary Ammonium | Quaternary Ammonium | Chlorine Dioxide |
| FDA Approval | | N | N | N | N | N | Y | N |
| Disinfectant Type | | Wipes | Liquid | Wipes | Wipes | Wipes | Spray | Spray |
| Convex | C1-6/C1-6T/C1-6i/SC1-6 /SC1-4HS/SC1-4H/ SC1-6H/ SC2-9H/SC1-7H | Δ | Δ | O | O | O | O | O |
| | C1-6CT | Δ | Δ | Δ | Δ | O | Δ | O |
| Micro-Convex | C5-8 | O | O | Δ | Δ | Δ | Δ | Δ |
| | C5-8N/C5-8NT/SC2-11H | Δ | Δ | O | O | O | O | O |
| Linear | L3-8/L3-8i/L3-12/L3-12T /L3-12i/L3-12H/ L3-12H ^{WD} /L3-2X/L8-17 /L8-17i/L8-17X/L3-15H | Δ | Δ | O | O | O | O | O |
| | L10-25H/SL3-19H /SL3-19X | Δ | Δ | O | O | O | Δ | O |
| Phased | MP1-5X/ SP1-5/SP1-5i/ SP1-5X/P1-5CT/SP1-5T/ SP3-8/SP3-8T | Δ | Δ | O | O | O | O | O |
| Volume | VC1-6/VC1-6T/SVC1-6/ SVC1-6H/VE3-10H/ SVC1-8H | Δ | Δ | O | O | O | Δ | O |
| Endo | EV3-10/EV3-10T/ EV3-10H/EV3-10X | Δ | Δ | O | O | O | O | O |
| | EC3-10/EC3-10T/ EC3-10H/EC3-10X/ EC(V)2-11H | Δ | Δ | O | O | O | O | O |
| | EN3-10 | Δ | Δ | O | Δ | Δ | Δ | Δ |
| | E3-10 | Δ | Δ | O | O | O | Δ | Δ |
| | E3-10H | Δ | Δ | O | O | O | Δ | Δ |
| IO | IO3-12/IO8-17/IO8-17T | Δ | Δ | O | O | O | O | O |
| CW | CW2.0/CW5.0 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| TEE | TEE3-7 | Δ | O | Δ | Δ | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | Δ | Δ | O | O | O | O | O |
| Others | US-Guider | Δ | Δ | O | O | O | O | X |
| | Strain relief&Cable | Δ | Δ | O | O | O | X | O |

※ CAUTION ※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

***NOTE**

- Antigermix, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant

O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Disinfectant | | | | | | |
|--------------------|--|---------------------------------|--------------------|-----------------------|----------------------------|-------------|-------------|-----------------------|
| | | SteriMax® Sporicide Wipes | SteriMax® Wipes | SteriMax Sporicide | UltraSan Broad Spectrum | CaviCide | CaviWipes | Metricide OPA Plus |
| Manufacturer | | Aseptix | Aseptix | Aseptix | Aseptix | Metrex | Metrex | Metrex |
| Level | | - | - | - | - | - | LLD,ILD | HLD |
| Active Ingredients | | Hydrogen peroxide | Hydrogen peroxide | Hydrogen peroxide | Hydrogen peroxide | Isopropanol | Isopropanol | Orthophthalaldehyde |
| FDA Approval | | N | N | N | N | Y | Y | Y |
| Disinfectant Type | | Wipes | Wipes | Spray | Spray | Spray | Wipes | Liquid |
| Convex | C1-6/C1-6T/C1-6i/SC1-6/SC1-4HS/SC1-4H/SC1-6H/ SC2-9H/SC1-7H | O | X | X | X | X | O | O |
| | C1-6CT | O | O | Δ | Δ | O | O | Δ |
| Micro-Convex | C5-8 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| | C5-8N/C5-8NT/SC2-11H | O | X | X | X | X | O | O |
| Linear | L3-8/L3-8i/L3-12/L3-12T/L3-12i/L3-12H/L3-12H ^{WD} /L3-2X/L8-17/L8-17i/L8-17X/L3-15H | O | X | X | X | X | O | O |
| | L10-25H/SL3-19H/SL3-19X | Δ | Δ | Δ | Δ | Δ | Δ | O |
| Phased | MP1-5X/ SP1-5/SP1-5i/SP1-5X/P1-5CT/SP1-5T/SP3-8/SP3-8T | O | X | X | X | X | O | O |
| Volume | VC1-6/VC1-6T/SVC1-6/SVC1-6H/VE3-10H/SVC1-8H | Δ | Δ | Δ | Δ | Δ | Δ | O |
| Endo | EV3-10/EV3-10T/EV3-10H/EV3-10X | O | X | X | X | X | O | O |
| | EC3-10/EC3-10T/EC3-10H/EC3-10X/EC(V)2-11H | O | X | X | X | X | O | O |
| | EN3-10 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| | E3-10 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| | E3-10H | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| IO | IO3-12/IO8-17/IO8-17T | O | X | X | X | X | O | O |
| CW | CW2.0/CW5.0 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| TEE | TEE3-7 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | O | X | X | X | X | O | O |
| Others | US-Guider | X | X | X | X | X | O | O |
| | Strain relief&Cable | O | O | X | X | O | O | O |

※ CAUTION ※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

***NOTE**

- Antigermix, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant

O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Disinfectant | | | | | | |
|--------------------|---|----------------|----------------------------|------------------------|-------------------|----------------------|-------------------|-------------------|
| | | Metricide | Meliseptol Wipes sensitive | Meliseptol Wipes ultra | Revital-OX Resert | Meliseptol Foam pure | Sporox II | Protex wipes |
| Manufacturer | | Metrex | B. Braun | B. Braun | STERIS | B. Braun | SultanHealthcare | PARKER |
| Level | | HLD | - | - | HLD | - | HLD | - |
| Active Ingredients | | Glutaraldehyde | Propan-1-ol | Quaternary Ammonium | Hydrogen Peroxide | Propan-1-ol | Hydrogen Peroxide | Ammonium chloride |
| FDA Approval | | Y | N | N | Y | N | Y | Y |
| Disinfectant Type | | Liquid | Wipes | Wipes | Liquid | Liquid | Liquid | Wipes |
| Convex | C1-6/C1-6T/C1-6i/SC1-6 /SC1-4HS/SC1-4H/ SC1-6H/ SC2-9H/SC1-7H | O | X | O | O | O | O | O |
| | C1-6CT | Δ | Δ | Δ | O | O | X | O |
| Micro-Convex | C5-8 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| | C5-8N/C5-8NT/SC2-11H | O | X | O | O | O | O | O |
| Linear | L3-8/L3-8i/L3-12/L3-12T /L3-12i/L3-12H/ L3-12H ^{WD} /L3-2X/L8-17 /L8-17i/L8-17X/L3-15H | O | X | O | O | O | O | O |
| | L10-25H/SL3-19H /SL3-19X | Δ | O | O | O | O | O | O |
| Phased | MP1-5X/ SP1-5/SP1-5i/ SP1-5X/P1-5CT/SP1-5T/ SP3-8/SP3-8T | O | X | O | O | O | O | O |
| Volume | VC1-6/VC1-6T/SVC1-6/ SVC1-6H/VE3-10H/ SVC1-8H | Δ | O | O | O | O | O | O |
| Endo | EV3-10/EV3-10T/ EV3-10H/EV3-10X | O | X | O | O | O | O | O |
| | EC3-10/EC3-10T/ EC3-10H/EC3-10X/ EC(V)2-11H | O | X | O | O | O | O | O |
| | EN3-10 | Δ | Δ | Δ | Δ | Δ | X | Δ |
| | E3-10 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| | E3-10H | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| IO | IO3-12/IO8-17/IO8-17T | O | X | O | O | O | O | O |
| CW | CW2.0/CW5.0 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| TEE | TEE3-7 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | O | X | O | O | O | O | O |
| Others | US-Guider | X | Δ | Δ | Δ | Δ | Δ | Δ |
| | Strain relief&Cable | O | O | O | X | O | Δ | Δ |

※ CAUTION ※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

***NOTE**

- Antigermix, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant

O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Disinfectant | | | | | |
|--------------------|--|---------------------|-------------------------------|-----------------------------|------------------|------------------|----------------------------|
| | | Sani-cloth AF3 | Acrylan | Descogen Liquid | Cutasept F | Cutasept G | Mikrobac Virucidal Tissues |
| Manufacturer | | NICE PAK | ANTISEPTICA | ANTISEPTICA | BODE Chemie GmbH | BODE Chemie GmbH | BODE Chemie GmbH |
| Level | | LLD,ILD | - | - | - | - | - |
| Active Ingredients | | Quaternary Ammonium | Quaternary ammonium compounds | Potassium peroxymonosulfate | Propan-2-ol | Propan-2-ol | ammonium chloride |
| FDA Approval | | Y | N | N | N | N | N |
| Disinfectant Type | | Wipes | Liquid | Liquid | Liquid | Liquid | Wipes |
| Convex | C1-6/C1-6T/C1-6i/SC1-6/SC1-4HS/SC1-4H/SC1-6H/SC2-9H/SC1-7H | O | X | O | X | X | O |
| | C1-6CT | O | Δ | Δ | X | X | O |
| Micro-Convex | C5-8 | Δ | Δ | Δ | Δ | Δ | Δ |
| | C5-8N/C5-8NT/SC2-11H | O | X | O | X | X | O |
| Linear | L3-8/L3-8i/L3-12/L3-12T/L3-12i/L3-12H/L3-12H ^{WD} /L3-2X/L8-17/L8-17i/L8-17X/L3-15H | O | X | O | X | X | O |
| | L10-25H/SL3-19H/SL3-19X | O | O | O | O | Δ | O |
| Phased | MP1-5X/ SP1-5/SP1-5i/SP1-5X/P1-5CT/SP1-5T/SP3-8/SP3-8T | O | X | O | X | X | O |
| Volume | VC1-6/VC1-6T/SVC1-6/SVC1-6H/VE3-10H/SVC1-8H | O | O | O | O | Δ | O |
| Endo | EV3-10/EV3-10T/EV3-10H/EV3-10X | O | X | O | X | X | O |
| | EC3-10/EC3-10T/EC3-10H/EC3-10X/EC(V)2-11H | O | X | O | X | X | O |
| | EN3-10 | Δ | Δ | Δ | Δ | Δ | Δ |
| | E3-10 | Δ | Δ | Δ | Δ | Δ | Δ |
| | E3-10H | Δ | Δ | Δ | Δ | Δ | Δ |
| IO | IO3-12/IO8-17/IO8-17T | O | X | O | X | X | O |
| CW | CW2.0/CW5.0 | O | Δ | Δ | Δ | Δ | Δ |
| TEE | TEE3-7 | Δ | Δ | Δ | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | O | X | O | X | X | O |
| Others | US-Guider | Δ | Δ | Δ | Δ | Δ | Δ |
| | Strain relief&Cable | Δ | Δ | Δ | Δ | Δ | Δ |

※CAUTION※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

*NOTE

- Antigermix, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant

O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Disinfectant | | | | |
|--------------------|--|---------------------------------------|----------------------------|-------------------------------|-------------------|-------------------|
| | | SONO Ultrasound and Disinfecting Wipe | Bacillol 30 Sensitive Foam | Bacillol 30 Sensitive Tissues | Incidin OxyFoam S | Incidin OxyWipe S |
| Manufacturer | | ADVANCED ULTRASOUND SOLUTIONS INC. | BODE Chemie GmbH | BODE Chemie GmbH | ECOLAB | ECOLAB |
| Level | | - | - | - | LLD,ILD | LLD,ILD |
| Active Ingredients | | ammonium chloride | Ethanol Propan-2-ol | Ethanol Propan-2-ol | Hydrogen peroxide | Hydrogen peroxide |
| FDA Approval | | N | N | N | N | N |
| Disinfectant Type | | Wipes | Liquid | Wipes | Liquid | Wipes |
| Convex | C1-6/C1-6T/C1-6i/SC1-6/SC1-4HS/SC1-4H/SC1-6H/ SC2-9H/SC1-7H | X | O | O | X | O |
| | C1-6CT | O | O | O | X | O |
| Micro-Convex | C5-8 | Δ | Δ | Δ | Δ | Δ |
| | C5-8N/C5-8NT/SC2-11H | X | O | O | X | O |
| Linear | L3-8/L3-8i/L3-12/L3-12T/L3-12i/L3-12H/L3-12H ^{WD} /L3-2X/L8-17/L8-17i/L8-17X/L3-15H | X | O | O | X | O |
| | L10-25H/SL3-19H/SL3-19X | O | O | O | O | O |
| Phased | MP1-5X/ SP1-5/SP1-5i/SP1-5X/P1-5CT/SP1-5T/SP3-8/SP3-8T | X | O | O | X | O |
| Volume | VC1-6/VC1-6T/SVC1-6/SVC1-6H/VE3-10H/SVC1-8H | O | O | O | O | O |
| Endo | EV3-10/EV3-10T/EV3-10H/EV3-10X | X | O | O | X | O |
| | EC3-10/EC3-10T/EC3-10H/EC3-10X/EC(V)2-11H | X | O | O | X | O |
| | EN3-10 | Δ | Δ | Δ | Δ | Δ |
| | E3-10 | Δ | Δ | Δ | Δ | Δ |
| | E3-10H | Δ | Δ | Δ | Δ | Δ |
| IO | IO3-12/IO8-17/IO8-17T | X | O | O | X | O |
| CW | CW2.0/CW5.0 | Δ | Δ | Δ | Δ | Δ |
| TEE | TEE3-7 | Δ | Δ | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | X | O | O | X | O |
| Others | US-Guider | Δ | Δ | Δ | Δ | Δ |
| | Strain relief&Cable | Δ | Δ | Δ | X | O |

※CAUTION※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

***NOTE**

- Antigermix, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant

O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Disinfectant | | | | |
|--------------------|---|--|-----------------------------------|---------------------|---------------------------------------|---|
| | | mikrozid® sensitive wipes | mikrozid® universal wipes premium | mikrozid® PAA wipes | gigasept® PAA concentrate | gigasept® pearls |
| Manufacturer | | schuelke | schuelke | schuelke | schuelke | schuelke |
| Level | | LLD, ILD | LLD, ILD | LLD, ILD | HLD | HLD |
| Active Ingredients | | Alkyl dimethylbenzyl ammonium chloride Didecylmethylammonium chloride Alkyl ethylbenzylammonium chloride | Propan-2-ol Ethanol | peracetic acid | peracetic acid potassium hydroxide | Sodium percarbonate Tetraacetylenedia mine |
| FDA Approval | | N | N | N | N | N |
| Disinfectant Type | | Wipes | Wipes | Wipes | Liquid | Powder |
| Convex | C1-6/C1-6T/C1-6i/SC1-6 /SC1-4HS/SC1-4H/ SC1-6H/ SC2-9H/SC1-7H | O | O | X | O | O |
| | C1-6CT | Δ | Δ | Δ | Δ | Δ |
| Micro-Convex | C5-8 | Δ | Δ | Δ | Δ | Δ |
| | C5-8N/C5-8NT/SC2-11H | O | O | X | O | O |
| Linear | L3-8/L3-8i/L3-12/L3-12T /L3-12i/L3-12H/ L3-12H ^{WD} /L3-2X/L8-17 /L8-17i/L8-17X/L3-15H | O | O | X | O | O |
| | L10-25H/SL3-19H /SL3-19X | O | O | X | O | O |
| Phased | MP1-5X/ SP1-5/SP1-5i/ SP1-5X/P1-5CT/SP1-5T/ SP3-8/SP3-8T | O | O | X | O | O |
| Volume | VC1-6/VC1-6T/SVC1-6/ SVC1-6H/VE3-10H/ SVC1-8H | O | O | X | O | O |
| Endo | EV3-10/EV3-10T/ EV3-10H/EV3-10X | O | O | X | O | O |
| | EC3-10/EC3-10T/ EC3-10H/EC3-10X/ EC(V)2-11H | O | O | X | O | O |
| | EN3-10 | Δ | Δ | Δ | Δ | Δ |
| | E3-10 | Δ | Δ | Δ | Δ | Δ |
| | E3-10H | Δ | Δ | Δ | Δ | Δ |
| IO | IO3-12/IO8-17/IO8-17T | O | O | X | O | O |
| CW | CW2.0/CW5.0 | Δ | Δ | Δ | Δ | Δ |
| TEE | TEE3-7 | Δ | Δ | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | O | O | X | O | O |
| Others | US-Guider | Δ | Δ | Δ | Δ | Δ |
| | Strain relief&Cable | O | O | X | O | O |

※CAUTION※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

***NOTE**

- Antigermix, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Int HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant
 O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Pre-Cleaner | | | | | |
|--------------------|--|---|---|----------------------------------|---------------------|-------------------|------------------|
| | | T-spray | T-spray II | Enzol | Theracide Plus(S5N) | Transeptic | Metrizyme |
| Manufacturer | | Pharmaceutical Innovations, INC | Pharmaceutical Innovations, INC | Johnson and Johnson | Burnishine Products | PARKER | Metrex |
| Level | | - | - | - | - | - | - |
| Active Ingredients | | Alkyl dimethyl benzyl ammonium chloride | Alkyl dimethyl benzyl ammonium chloride | Subtilisins(proteolytic enzymes) | Isopropyl alcohol | Isopropyl alcohol | Propylene glycol |
| FDA Approval | | Y | Y | Y | N | Y | Y |
| Disinfectant Type | | Spray | Spray | Liquid | Liquid | Liquid | Liquid |
| Convex | C1-6/C1-6T/C1-6i/SC1-6/SC1-4HS/SC1-4H/SC1-6H/ SC2-9H/SC1-7H | Δ | Δ | O | Δ | Δ | O |
| | C1-6CT | Δ | Δ | Δ | Δ | Δ | O |
| Micro-Convex | C5-8 | O | O | Δ | Δ | Δ | Δ |
| | C5-8N/C5-8NT/SC2-11H | Δ | Δ | O | Δ | Δ | O |
| Linear | L3-8/L3-8i/L3-12/L3-12T/L3-12i/L3-12H/L3-12H ^{WD} /L3-2X/L8-17/L8-17i/L8-17X/L3-15H | Δ | Δ | O | Δ | Δ | O |
| | L10-25H/SL3-19H/SL3-19X | Δ | O | O | O | O | O |
| Phased | MP1-5X/ SP1-5/SP1-5i/SP1-5X/P1-5CT/SP1-5T/SP3-8/SP3-8T | Δ | Δ | O | Δ | Δ | O |
| Volume | VC1-6/VC1-6T/SVC1-6/SVC1-6H/VE3-10H/SVC1-8H | Δ | O | O | O | O | O |
| Endo | EV3-10/EV3-10T/EV3-10H/EV3-10X | Δ | Δ | O | Δ | Δ | O |
| | EC3-10/EC3-10T/EC3-10H/EC3-10X/EC(V)2-11H | Δ | Δ | O | Δ | Δ | O |
| | EN3-10 | O | Δ | Δ | Δ | Δ | O |
| | E3-10 | O | Δ | Δ | Δ | Δ | Δ |
| | E3-10H | O | Δ | Δ | Δ | Δ | Δ |
| IO | IO3-12/IO8-17/IO8-17T | Δ | Δ | O | Δ | Δ | O |
| CW | CW2.0/CW5.0 | Δ | Δ | Δ | Δ | Δ | Δ |
| TEE | TEE3-7 | Δ | Δ | O | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | Δ | Δ | O | Δ | Δ | O |
| Others | US-Guider | Δ | Δ | O | Δ | O | X |
| | Strain relief&Cable | Δ | Δ | O | Δ | O | O |

※ CAUTION ※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

***NOTE**

- Antigermix, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant

O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Pre-Cleaner | | | | | |
|--------------------|--|-------------|----------------------|---------------------|-------------------|------------------|---------------------------------|
| | | Klenzyme | McKesson | Cidezyme | Endozyme | EmPower | Ultrasound probe cleaning wipes |
| Manufacturer | | STERIS | McKesson | Johnson and Johnson | RUHOF Corporation | Metrex | Nanosonics |
| Level | | - | - | - | - | - | - |
| Active Ingredients | | Enzymes | PCMX (Chloroxylenol) | Proteolytic Enzymes | Propanol | Propylene glycol | PHENOXYETHANOL |
| FDA Approval | | Y | Y | Y | Y | Y | N |
| Disinfectant Type | | Liquid | Liquid | Liquid | Liquid | Liquid | Wipes |
| Convex | C1-6/C1-6T/C1-6i/SC1-6/SC1-4HS/SC1-4H/SC1-6H/ SC2-9H/SC1-7H | Δ | Δ | Δ | Δ | O | O |
| | C1-6CT | Δ | Δ | Δ | Δ | O | Δ |
| Micro-Convex | C5-8 | Δ | Δ | Δ | Δ | Δ | Δ |
| | C5-8N/C5-8NT/SC2-11H | Δ | Δ | Δ | Δ | O | O |
| Linear | L3-8/L3-8i/L3-12/L3-12T/L3-12i/L3-12H/L3-12H ^{WD} /L3-2X/L8-17/L8-17i/L8-17X/L3-15H | Δ | Δ | Δ | Δ | O | O |
| | L10-25H/SL3-19H/SL3-19X | Δ | Δ | Δ | O | O | O |
| Phased | MP1-5X/ SP1-5/SP1-5i/SP1-5X/P1-5CT/SP1-5T/SP3-8/SP3-8T | Δ | Δ | Δ | Δ | O | O |
| Volume | VC1-6/VC1-6T/SVC1-6/SVC1-6H/VE3-10H/SVC1-8H | Δ | Δ | Δ | O | O | O |
| Endo | EV3-10/EV3-10T/EV3-10H/EV3-10X | Δ | Δ | Δ | Δ | O | O |
| | EC3-10/EC3-10T/EC3-10H/EC3-10X/EC(V)2-11H | Δ | Δ | Δ | Δ | O | O |
| | EN3-10 | O | O | O | Δ | Δ | Δ |
| | E3-10 | Δ | Δ | Δ | Δ | Δ | Δ |
| | E3-10H | Δ | Δ | Δ | Δ | Δ | Δ |
| IO | IO3-12/IO8-17/IO8-17T | Δ | Δ | Δ | Δ | O | O |
| CW | CW2.0/CW5.0 | Δ | Δ | Δ | Δ | Δ | Δ |
| TEE | TEE3-7 | Δ | Δ | O | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | Δ | Δ | Δ | Δ | O | O |
| Others | US-Guider | Δ | Δ | Δ | Δ | O | Δ |
| | Strain relief&Cable | Δ | Δ | Δ | Δ | O | Δ |

※ CAUTION ※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

***NOTE**

- Antigermin, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant

O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Gel | | | | | Trophon® EPR | |
|--------------------|--|--------------|----------------------------|-----------------------------|--------|------------------------|-----------------|------------------------|
| | | Aquasonic100 | ULTRASOUND GEL | Ultrasound Transmission gel | Scan | Sonogel | | Super sonic |
| Manufacturer | | PARKER | Ultragel Hungary 2000 Ltd. | Ecosonic | PARKER | SONOGEL Vertriebs GmbH | | Sung Heung Corporation |
| Level | | - | - | - | - | - | | - |
| Active Ingredients | | - | - | - | - | - | | - |
| FDA Approval | | N | N | N | N | N | | N |
| Disinfectant Type | | Gel | Gel | Gel | Gel | Gel | | Gel |
| Device | | | | | | | | |
| Convex | C1-6/C1-6T/C1-6i/SC1-6/SC1-4HS/SC1-4H/SC1-6H/ SC2-9H/SC1-7H | O | O | O | Δ | O | O | O |
| | C1-6CT | O | Δ | Δ | Δ | Δ | Δ | Δ |
| Micro-Convex | C5-8 | O | O | Δ | O | O | Δ | Δ |
| | C5-8N/C5-8NT/SC2-11H | O | O | O | Δ | O | O | O |
| Linear | L3-8/L3-8i/L3-12/L3-12T/L3-12i/L3-12H/L3-12H ^{WD} /L3-2X/L8-17/L8-17i/L8-17X/L3-15H | O | O | O | Δ | O | O | O |
| | L10-25H/SL3-19H/SL3-19X | O | Δ | Δ | Δ | O | Δ | O |
| Phased | MP1-5X/ SP1-5/SP1-5i/ SP1-5X/P1-5CT/SP1-5T/ SP3-8/SP3-8T | O | O | O | Δ | O | O | O |
| Volume | VC1-6/VC1-6T/SVC1-6/ SVC1-6H/VE3-10H/ SVC1-8H | O | Δ | Δ | Δ | O | Δ | O |
| Endo | EV3-10/EV3-10T/ EV3-10H/EV3-10X | O | O | O | Δ | O | O | X |
| | EC3-10/EC3-10T/ EC3-10H/EC3-10X/ EC(V)2-11H | O | O | O | Δ | O | O | O |
| | EN3-10 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| | E3-10 | O | O | Δ | Δ | Δ | Δ | Δ |
| | E3-10H | O | O | Δ | Δ | Δ | Δ | Δ |
| IO | IO3-12/IO8-17/IO8-17T | O | O | O | Δ | O | O | O |
| CW | CW2.0/CW5.0 | Δ | Δ | Δ | Δ | Δ | Δ | Δ |
| TEE | TEE3-7 | O | Δ | Δ | Δ | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | O | O | O | Δ | O | O | X |
| Others | US-Guider | O | Δ | Δ | Δ | Δ | X | Δ |
| | Strain relief&Cable | O | Δ | Δ | Δ | O | O | Δ |

※CAUTION※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

***NOTE**

- Antigermix, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25

Transducer Disinfectant Material Compatibility Table

[Rev.25]

HLD : High level disinfectant / ILD : Intermediate level disinfectant / LLD : Low level disinfectant

O : Approved / X : Not Approved / Δ : Not Tested

| Product Type | | Trophon 2 | Antigermix ^{Note)} | V-PRO ^{Note)} |
|--------------------|---|-----------|-----------------------------|------------------------|
| Manufacturer | | | | |
| Level | | | | |
| Active Ingredients | | | | |
| FDA Approval | | | | |
| Disinfectant Type | | Device | Device | Device |
| Convex | C1-6/C1-6T/C1-6i/SC1-6 /SC1-4HS/SC1-4H/ SC1-6H/ SC2-9H/SC1-7H | O | O | O |
| | C1-6CT | Δ | Δ | O |
| Micro-Convex | C5-8 | Δ | Δ | Δ |
| | C5-8N/C5-8NT/SC2-11H | O | O | O |
| Linear | L3-8/L3-8i/L3-12/L3-12T /L3-12i/L3-12H/ L3-12H ^{WD} /L3-2X/L8-17 /L8-17i/L8-17X/L3-15H | O | O | O |
| | L10-25H/SL3-19H /SL3-19X | O | X | O |
| Phased | MP1-5X/ SP1-5/SP1-5i/ SP1-5X/P1-5CT/SP1-5T/ SP3-8/SP3-8T | O | O | O |
| Volume | VC1-6/VC1-6T/SVC1-6/ SVC1-6H/VE3-10H/ SVC1-8H | O | X | O |
| Endo | EV3-10/EV3-10T/ EV3-10H/EV3-10X | O | O | O |
| | EC3-10/EC3-10T/ EC3-10H/EC3-10X/ EC(V)2-11H | O | O | O |
| | EN3-10 | Δ | Δ | Δ |
| | E3-10 | Δ | Δ | Δ |
| | E3-10H | Δ | Δ | Δ |
| IO | IO3-12/IO8-17/IO8-17T | O | O | O |
| CW | CW2.0/CW5.0 | Δ | Δ | Δ |
| TEE | TEE3-7 | Δ | Δ | Δ |
| minisono | minisono L3-12/C1-6 | X | O | X |
| Others | US-Guider | Δ | Δ | Δ |
| | Strain relief&Cable | Δ | X | O |

※CAUTION※

- Only wipe the cable, strain relief and parts of the transducer. Do not wipe the connector of the transducer.
- Do not immerse transducers in a solution for longer than one hour, unless under disinfectant manufacturer's instructions.
- Using a non-recommended disinfectant or not following the recommended disinfection method can damage and/or discolor the transducer and will void the transducer warranty.

***NOTE**

- Antigermix, V-PRO : Slight discoloration of the housing and cable can occur. These changes have no impact on device performance

Alpinion Medical Systems Co., Ltd.

e-mail : marketing@alpinion.com TEL : +82-2-3282-0900

Copyright© 2017-2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved

MKT-EL-TS-LF-E2207-V25